

BIO365L TA Cleanup Checklist

Date: _____

Check off Student Cleanup

- ____ Equipment turned off on all rigs
- ____ Vacuum flasks empty
- ____ Computers logged off

Lab Shutdown

- ____ Rinse slice holding gas diffuser
- ____ Turn oxygen cylinders off
- ____ Turn vacuum pump off
- ____ Empty slice holding beaker & rinse
- ____ Vibratome tray cleaned & rinsed
- ____ Vibratome specimen disc cleaned & rinsed
- ____ Cutting saline beaker rinsed & dried
- ____ Dissection tools cleaned & stored
- ____ Animal carcass to freezer
- ____ Erase whiteboard
- ____ Check stock of fire-polished pipette glass
- ____ Check stock of pipette filler tips